# TRUEread® Blood Glucose Test Strips

## Instructions for Use

### Intended Use

TRUEread Blood Glucose Test Strips are used with TRUEread or TrackEASE Meters to quantitatively measure whole blood glucose. The TRUEread can be used for self-testing at home or for professional use. The TRUEread is for *in vitro* (external) use only.

### Test Principle

TRUEread measures glucose by using amperometric technology employing a glucose oxidase reaction. When whole blood or Control is drawn into the tip of a Test Strip, glucose in the sample reacts with chemicals and produces an electrical current. The Meter measures electrical current and calculates amount of glucose. The glucose result is displayed as a plasma value.

### Chemical Composition

Glucose oxidase (*Aspergillus* sp.) 2.5 units, mediators, buffers and stabilizers.

① **Contact End** End of Strip inserted into Meter.

② **Sample Tip** Pointed end of Strip where sample is drawn into Strip.

Top of Test Strip

②

Correct

Incorrect

### Caring for Test Strips

- Use Test Strips quickly after removing from vial. Recap vial right away. Test Strips left outside of vial too long may give an error message
- Test Strips **must** be kept in original vial, with vial cap tightly sealed. NEVER transfer Strips from one vial to another.
- Write date opened on vial label when removing first Strip. Discard all unused Strips in vial after date printed on the Strip vial label, or 120 days after date opened, if either date has passed.
- Store Test Strip vials in a dry place at room temperature below 86°F/30°C.  
**DO NOT REFRIGERATE OR FREEZE.**
- DO NOT store in bathroom or kitchen.
- DO NOT expose to extreme heat or cold, direct sunlight, or high humidity for any length of time.
- DO NOT reuse Test Strips.
- DO NOT bend, cut or alter Test Strips in any way.

### Important Information

- Test Strips are for *in vitro* testing only. **DO NOT consume.**
- ONLY** use TRUEread Test Strips with TRUEread or TrackEASE Meters and TRUEcontrol Glucose Control. Using other Meters or Control may cause inaccurate results.
- Do not change your medication or treatment plan based on results from the TRUEread Strips and TRUEread or TrackEASE Meters without advice of Doctor or Diabetes Healthcare Professional.
- Discard used Test Strips and Lancets into an appropriate container (biohazard/for sharps objects).
- Using Test Strips past written or printed dates on the Strip vial label may cause inaccurate results.
- Discard any strips or vials that appear damaged.
- NEVER** use serum, plasma, or clotted blood when testing with TRUEread. **ONLY** use fresh capillary whole blood from the finger or forearm.
- TRUEread is not recommended for use with venous samples.**
- When using the forearm for blood sample:<sup>1</sup>
  - Check with your Doctor or Diabetes Healthcare Professional to see if forearm testing is right for you.
  - Results from forearm are not always the same as results from finger.
- Use finger instead of forearm for more accurate results:
  - Within 2 hours of eating, exercise, or taking insulin,
  - If your blood sugar may be rising or falling rapidly or your routine results are often fluctuating,
  - If you are ill or under stress,
  - If your forearm test results do not match how you feel,
  - If your blood sugar may be low or high,
  - If you do not notice symptoms when blood sugar is high or low.

### Quality Control (QC) Testing

There are three quality control tests to let you know that the System is working properly. An automatic self-test is performed each time a Strip is inserted into a Meter. Upon inserting a Strip into the Test Port, if all segments appear and the Apply Sample Symbol appears in the Display, the Meter is working properly.

### Quality Control Test: Check Strip

TRUEread Check Strip is used to make sure that the Meter is working properly. The Check Strip is a hard plastic strip found in Carrying Case of TRUEread Meter. When testing with Check Strip, OK~ is shown in Meter Display if Meter is working properly. See Owner's Booklet for more information on Check Strip testing.

### Quality Control Test: Control

TRUEcontrol Glucose Control is used to check testing technique and System performance. When Control results fall within ranges found on Test Strip vial label of Test Strips being used, System is working properly and testing technique is good. See Owner's Booklet or TRUEcontrol Glucose Control Instructions for Use for more information on Control Testing.

### Important Information

There are three levels of TRUEcontrol Glucose Control available that contain known amounts of glucose. It is important to perform Control Tests with more than one level of Control to assure your System is working properly and your testing technique is good. For more information on obtaining different levels of control, call 1-800-803-6025 or 1-954-677-4599.

### Blood Glucose Testing

- Wash area to be lanced. Dry thoroughly.
- Allow Meter and Test Strips to sit at room temperature for 30 minutes before using. If opening vial for the first time, write date opened on vial label.
- Check EXP Date. Do not use if either date has passed. Discard vial and test with new vial.
- Remove one Strip from vial. Recap vial immediately.
- Insert Contact End of Strip, contacts facing up, into Test Port of Meter. Meter turns on. Do not remove Strip from Meter until testing is finished.
- Match Code in Meter Display to Code on Strip vial being used. If Codes do not match, see *Troubleshooting*. Do not use Strips if Codes do not match.
- Obtain blood drop.
- After Apply Sample Symbol appears in Display place Sample Tip of Strip gently against blood drop and hold. Allow Strip to draw in blood (1 µl) until meter beeps. Remove Strip from drop. Dashes appear as Meter is testing. If Meter does not beep 5 seconds after touching Strip to blood drop, see *Troubleshooting*.
- Result is displayed. Record result.
- Remove Strip from Meter. Meter turns off. Discard used Strip and lancet in appropriate container.

### Expected Results for people without diabetes:<sup>2</sup>

	Plasma Blood Glucose Result
Before eating	<100 mg/dL
Two hours after meals	<140 mg/dL
Bedtime	<120 mg/dL

A Doctor or Diabetes Healthcare Professional determines personal target glucose ranges. Always check with your Doctor or Diabetes Healthcare Professional before changing your treatment program. If you are having symptoms that suggest your glucose is too low or too high, contact your Doctor or Diabetes Healthcare Professional right away. The TRUEread is optimized for testing with fresh capillary whole blood and calibrated to display a plasma result. Plasma results are estimated to be 10%-15% higher than whole blood results.<sup>3</sup> If comparing results using TRUEread Strips to laboratory test results, perform a fingerstick blood test within 15 minutes of the laboratory test. Results are considered accurate when they are within 20% of laboratory results. If you have eaten recently, results using TRUEread Strips can be up to 70 mg/dL higher than laboratory results. If glucose tests are performed using capillary blood with the TRUEread and venous blood with the laboratory system, capillary blood (finger) may be up to 70 mg/dL higher than venous blood (vein), even when each test is performed within 10 minutes after sampling and per manufacturers' instructions.<sup>3</sup>

### Troubleshooting (See Owner's Booklet *Display Messages* for more details.)

If your result is unusually high or low or doesn't match the way you feel, perform a Control Test (see *Quality Control Tests*).

If the Control Test is within range:

- Read Blood Glucose Testing again.
- Test your blood glucose again using a new TRUEread Strip.

If the results still do not match the way you feel, check with your Doctor or Diabetes Healthcare Professional before changing your treatment plan.

If the results are not within range:

- Check the Expiration Dates. Do not use if past either written date or date printed on Test Strip vial or Control bottle. Test with new Strips/Controls.
- Check for error messages. If an error message appears, follow the Actions in the Display Message Section of the Owner's Booklet.
- Check your testing technique. Perform another Control Test.
- Check the temperature. Allow System to reach room temperature 59°-86°F (15°-30°C) before testing.

### Limitations

TRUEread products cannot be used for testing blood glucose in newborns.

The following WILL affect accurate test results:<sup>4</sup>

- Abnormally high doses of acetaminophen.

The following MAY affect accurate test results:<sup>9,10</sup>

- Hyperglycemia with hyperosmolarity, with or without ketosis.
- Thickened blood caused by dehydration.

*NOTE:* Similar observations have been reported by other blood glucose systems.<sup>5,6</sup>

The following will NOT affect accurate test results:<sup>4</sup>

- Hematocrit levels between 30% and 55%.
- Elevated total blood cholesterol and triglycerides.
- Salicylate occurring in expected blood concentrations.
- Testing at altitudes up to and including 10,150 feet.

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patient should not be tested with blood glucose meters.

### Performance Characteristics<sup>2</sup>

#### Accuracy

Accuracy describes how well TRUEread results agree with results from an accepted laboratory analyzer. The accuracy of TRUEread was assessed against Yellow Springs Instrument Model 23A<sup>7</sup> (which also uses a glucose oxidase method). Testing was performed by lay users at 4 independent clinical sites. The study demonstrates that TRUEread compares well with the laboratory method. A regression analysis of the 229 results is given below:

N	-	229	Slope	-	0.98
Y-intercept	-	5.18 mg/dL	r	-	0.96

#### Precision

Precision describes the variation between TRUEread results. Precision studies were performed in a laboratory.

Within Run Precision (whole blood)

	<u>55 mg/dL</u>	<u>82 mg/dL</u>	<u>122 mg/dL</u>	<u>241 mg/dL</u>	<u>500 mg/dL</u>
C.V. (%)	5.1	4.2	3.7	2.7	3.0

Total Precision (Control Solution)

	<u>90 mg/dL</u>	<u>161 mg/dL</u>	<u>286 mg/dL</u>
C.V. (%)	6.3	5.6	4.5

#### Consumer Study

TRUEread was tested at diabetes clinics of a large urban area hospital and diabetes care centers. TRUEread results were obtained from 212 patients and were compared to parallel results obtained on Yellow Springs Model 2300 STAT. This study shows that patients, when self-testing, and laboratory professionals are able to obtain accurate results when compared to results obtained from a laboratory system.

A clinical study was performed at an urban area hospital to compare forearm and finger stick blood glucose results using the TRUEread. TRUEread blood glucose tests were performed by patients self-testing (n=100) and on patients' blood by healthcare professionals (n=108). Results demonstrate substantial equivalence between forearm and finger blood glucose values when patients' blood glucose is in steady state and blood glucose tests are performed by patients or healthcare professionals using the TRUEread.


**See TRUEread or TrackEASE Owner's Booklet for detailed instructions. Call Nipro Diagnostics, Inc. at 1-800-803-6025 (USA) or 1-954-677-4599 for assistance or visit our website at [www.niprodiagnostics.com](http://www.niprodiagnostics.com) for more information.**

**For medical assistance, contact a Doctor or Diabetes Healthcare Professional.**

#### References

- Food and Drug Administration 2006. *Diabetes Information - Glucose Meters and Diabetes Management* New Technologies: Alternative Site Testing [Electronic Version]. Retrieved February 17, 2006 from [www.fda.gov/diabetes/glucose.html#fll](http://www.fda.gov/diabetes/glucose.html#fll).
- Joslin Diabetes Center. *Goals for Blood Glucose Control* [Electronic Version]. Retrieved June 2, 2008 from [www.joslin.org/Beginners\\_guide\\_523.asp](http://www.joslin.org/Beginners_guide_523.asp)
- Sacks, D.B.; "Carbohydrates" Burtis, C.A. and Ashwood, E.R. (eds.), Tietz Textbook of Clinical Chemistry, Philadelphia: W.B. Saunders Company (1994), 959.
- Data on file.
- Young, D.S. Effects of Drugs of Clinical Laboratory Tests, Third Edition. Washington, D.C.: AACCC Press (1990), pp. 3-168, 3-180.
- Sacks, D.B., "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.
- Yellow Springs Instrument Model23A Glucose Analyzer Instruction Manual, "Principles of Operation", pp. 1-2 and "Performance Specifications", pp. 3, 1986.
- Atkins, S.H.; Dasamhapatra, A.; Jaker, M.A.; Choast, M.I.; Redd, S., Fingerstick Glucose Determination in Shock. *Annals of Internal Medicine*, 114:1020-1024 (1991).
- Sander, M.; Low-Ber, I., Misleading Capillary Glucose Measurements, *Practical Diabetes*, 7:210 (1990).
- Wickham, N.W.R.; Achor, K.N.; Cove, D.H., Unreliability of Capillary Blood Glucose in Peripheral Vascular Disease, *Practical Diabetes*, 3:100 (1986).

TRUEread

G3NP03R20

TRUEread

Manufactured by:

**NIPRO**  
LAUDNOSTICS..

Fort Lauderdale, FL 33309 U.S.A.  
[www.niprodiagnostics.com](http://www.niprodiagnostics.com)

TRUEread®

© 2010 Nipro Diagnostics, Inc. TrackEASE, TRUEread, TRUEcontrol, and the Nipro Diagnostics logo are trademarks of Nipro Diagnostics, Inc. G3NP03 Rev. 20

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