

TRUEbalance™ Blood Glucose Monitoring System

What is the TRUEbalance™ Blood Glucose Testing System?

Used to test blood glucose levels, TRUEbalance™ is a no coding meter that delivers outstanding performance at an affordable price; the TRUEbalance™ System is perfect for all patients.

What is no coding?

The TRUEbalance™ no coding system does not require you to perform an extra step to match up your meter with your test strips. Just insert a strip and test.

What strips do I use for the TRUEbalance™ meter?

TRUEbalance™ meters use TRUEbalance™ [Blood Glucose Test Strips](#).

How do I set up my TRUEbalance™ meter?

See [instructional video](#) and refer to your [owner's booklet](#) included with your meter.

How do I calibrate my TRUEbalance™ meter?

Coding is not necessary; TRUEbalance™ is a no code System.

How do I use my TRUEbalance™ System?

See [instructional video](#) and refer to your [owner's booklet](#) included with your meter.

Can the TRUEbalance™ System be used to test on alternative test sites?

Yes. TRUEbalance™ is suitable for fingertip and forearm testing.

Are special lancets required when testing with the TRUEbalance™ meter?

No. We recommend the TRUEdraw™ lancing device for greater comfort and control to be used with our single-use sterile lancets.

I am on peritoneal dialysis, can I use this meter?

Yes. The TRUEbalance™ System is a great choice for you. Please consult your instructions for use for limitations.

Does the TRUEbalance™ meter have memory? Does the meter store my results?

Yes. TRUEbalance™ has a 365-test memory with date and time and uploading capability.

Can I upload my meter data into data management software?

Yes. The TRUEbalance™ meter uses [TRUEmanager™ Diabetes Management Software](#). Follow the link to download the software or to request the software on CD. A SmartData Cable may be purchased separately.

Is the TRUEbalance™ meter covered by a lifetime warranty?

Yes. We are committed to ensuring that your meter and testing requirements are fulfilled. Our meters are covered by a lifetime warranty.

What kind of batteries does the TRUEbalance™ meter use?

The TRUEbalance™ meter uses a non rechargeable 3V lithium battery (#CR2032).

Where can I find the TRUEbalance™ meter?

Our TRUE portfolio of products is available at your favorite retail store or mail service provider. [Where to Buy](#)